

SHOPRIDER TE888SL

enjoy the freedom...

NEW UPGRADED MODEL
More Power • More Speed • More Life



This wonderfully comfortable scooter owes its super smooth ride to the 4 wheel suspension, large diameter wheels and its multi adjustable seat.

The big surprise however is its speed, powered by a 1.5 horsepower motor, you can now get to where you want to go in half the time.



Sign up to Shopriders AA Mobility Care programme and make sure you Get Home Safely with the AA

SHOPRIDER[®]
New Zealand

Authorised Dealer



FEATURES

- Powerful 1.5hp motor
- Headlights & reflectors for night use
- Modern ergonomically designed steering tiller
- On/Off security key
- Large range of accessories
- 4 wheel suspension

SPECIFICATIONS

width	600mm
length	1280mm
height (headrest lowered)	1100mm
weight (incl battery)	97kg
speed	0-11 kph
wheel size	279mm / 11"
max. climb rate	9° / 15%
turning radius	1460mm
max. load	136 kg

All specifications are correct as at April 2010 but are subject to change without notice. Maximum speed is subject to various conditions such as terrain, user weight and battery condition. Colour and options are subject to availability and may not be exactly as shown